



Accuris Mini Balances are extraordinarily compact, easy to use, and convenient for a wide range of laboratory weighing applications. Despite a footprint of only 2.75 x 4 inches, these Mini Balances offer exceptional accuracy and linearity through their full weighing ranges.

Both models feature an easy-to-read back lit liquid crystal display, and input buttons for ON/OFF, MODE, and TARE. The software includes a calibration function, and the balances can be easily and accurately calibrated using a 100g or 500g calibration weight.

The Mini Balances are powered with 2 x AAA batteries, included, and an auto shut off feature turns off the power after one minute of non-use to save battery life. Each balance includes a snap on, protective cover with instructions conveniently printed on the underside.

- Easy-to-read back lit liquid crystal display
- Calibration function
- Includes 2 x AAA batteries
- Snap on, protective cover



## Ordering Information & Technical Data:

Item No.	Description:
W4000-100	Accuris™ "Mini 100g" Balance, 100 grams/0.1g readability
W4000-500	Accuris™ "Mini 500g" Balance, 500 grams/0.1g readability

Model	W4000-100	W4000-500
Max Capacity	100 grams	500 grams
Readability	0.01 g	0.1 g
Repeatability	0.01 g	0.1 g
Linearity (+/-)	0.01 g	0.1 g
Calibration	External	
Pan Size	2.1 x 2.25 in / 5.3 x 5.7 cm	
Weighing units	g, ct, dwt, gn	
Display	backlit LCD	
Power	2 x AAA batteries	
Exterior Dims.	2.6 x 3.8 x 0.75 in / 6.5 x 9.7 x 1.8 cm	