

## Accuris W1000 Series Calibration Weights



Accuris™ calibration weights are ideal for calibration of your Accuris (or any brand) lab balance. Three single weights are available, 100g, 200g, and 500g, as well as a recently introduced complete weight set. All weights meet class F1 tolerance per OIML (International Organization of Legal Metrology) Recommendation R111 and are made of highly polished, non-magnetic stainless steel.

The W1000 series of individual weights come in protective plastic cases with a screw top cover. The W1100-200 Weight Set includes 23 weights ranging from 1mg to 200g and is packaged in a protective carrying case with a set of tweezers, a dust brush and cleaning cloth.



### Ordering Information:

W1000-100	Accuris™ Calibration Weight, 100g (F1 grade)
W1000-200	Accuris™ Calibration Weight, 200g (F1 grade)
W1000-500	Accuris™ Calibration Weight, 500g (F1 grade)
W1000-TF	Tweezers, for handling calibration weights
W1100-200	Calibration set, 23 pieces from 1mg to 200g with tweezers, dust brush, cleaning cloth & protective case

### Also available from Accuris Instruments:



Analytical, Precision and Compact Balances



Mini Lab Balances

Contact your Benchmark Scientific distributor to order.