

JoJo
Life Science



www.nestscientificusa.com

NEST[®]

2021

www.jojo-ls.de – info@jojo-ls.de

PCR Tubes & Caps


Features

- Free of Dnase and Rnase, human DNA and PCR inhibitors
- Autoclavable




PCR Tube

Cat.No.	Volume (mL)	/Pack	/Case
401001	0.2	1,000	10,000



PCR 8-strip Tubes, White

Cat.No.	Volume (mL)	/Pack	/Case
403112	0.1	125	1,250
403012	0.2	125	1,250




PCR 8-strip Tube Caps

Cat.No.	Cap Type	/Pack	/Case
406012	Flat	125	1,250




PCR 8-strip Tubes with Individual Caps Attached

Cat.No.	Volume (mL)	/Pack	/Case
404001	0.2	120	1,200



PCR 8-strip Tubes, Transparent

Cat.No.	Volume (mL)	/Pack	/Case
403102	0.1	125	1,250
403002	0.2	125	1,250



PCR 8-strip Tube Caps

Cat.No.	Cap Type	/Pack	/Case
406112	Domed	125	1,250



When using a thermal cycler with fluorescence detector above the samples, NEST white PCR plates effectively prevent crosstalk between wells. The Nest white PCR plate also allows for high accuracy of detection when signal intensity is low. Nest clear PCR plates are designed for PCR machines where the fluorescence detector is located below the sample plate.


Features

- Made of high quality polypropylene to ensure minimal loss of reaction solution.
- Flat surface, thick and solid, not easy to deform.
- Elevated edge of the hole can better prevent cross contamination.
- Compatible with NEST pressure sensitive film, self-adhesive film and hot sealing film.
- Autoclavable.

Compatible with Roche Systems



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402712	0.1	Semi Skirt	H12	10	50



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402812	0.1	No Skirt	A12/H12	10	50



Cat.No.	Volume (uL)	Skirt Type	Cut Position	/Pack	/Case
409013	40	Semi Skirt	A24/P24	10	50

PCR Adaptation Table

Brand	Instrument Model	Cat.No.		
		402712	402812	409013
Roche	lightcycler 480	✓		
	lightcycler 480 II	✓		✓
	lightcycler 96	✓		
	lightcycler nano		✓	

PCR Plate, No Skirt, Elevated Wells

Elevated wells plate is designed to reduce the risk of cross contamination.



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402201	0.2	No Skirt	H12	25	100

PCR Plate, Full Skirt

Stronger construction and more space for labeling.



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402501	0.1	Full Skirt	H1	25	100

PCR Plates, No Skirt


The low profile design minimizes the condensate contamination and evaporation of reaction solution.



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402101	0.1	No Skirt	H12	25	100



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402001	0.2	No Skirt	H1	25	100



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402111	0.1	No Skirt	H12	25	100



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402011	0.2	No Skirt	H1	25	100

PCR Plate, Semi Skirt, Compatible with ABI Machines

Stronger construction and more space for labeling.



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402401	0.1	Semi Skirt	A1	25	100



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402601	0.2	Semi Skirt	A12	25	100

PCR Plate, Semi Skirt

Stronger construction and more space for labeling.



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402301	0.2	Semi Skirt	A12	25	100



Cat.No.	Volume (mL)	Skirt Type	Cut Position	/Pack	/Case
402411	0.1	Semi Skirt	A1	25	100

PCR Rack

Features

- Fits all NEST PCR products
- Stackable
- Autoclavable

Cat.No.	Specification	/Pack	/Case
407001	8*12, One Rack & One Lid	1	25
407101	8*12, Five Racks & One Lid	5	25



PCR Sealing Film & Scraper



Cat.No.	Size (mm)	/Pack
411001	70*78	1



Cat.No.	Specification	/Pack	/Case
410001	Standard Sealing Film, -20°C To 110°C	100	500
410011	Advanced Sealing Film, -70°C To 100°C	100	500

Brand	Instrument Model	PCR Tube				PCR Plate							
		401001	403002/404001	401001	403002/404001	402001	402301	402401	402101	402601	402812		
Applied Biosystems®	2720	401001	403002/404001	401001	403002/404001	402001	402301				402601		
	9600	401001	403002/404001	401001	403002/404001	402001	402301				402601		
	9700	401001	403002/404001			402001	402301				402601		
	9800Fast				403102				402401	402101			
	Veriti/Fast				403102				402401	402101			
	7000	401001	403002/404001	401001	403002/404001	402001	402301				402601		
	7300	401001	403002/404001	401001	403002/404001	402001	402301				402601		
	7500/Fast	401001	403002/404001	401001	403002/404001	402001	402301		402401	402101		402601	
	7700	401001	403002/404001	401001	403002/404001	402001	402301					402601	
	7900HT/Fast	401001	403002/404001	401001	403002/404001	402001	402301		402401	402101		402601	
	StepOne				403102				402401	402101	402812		
	StepOnePlus				403102				402401	402101	402812		
	ProflexPCRSytem	401001	403002/404001	401001	403002/404001	402001	402301					402601	
	QuantStudio 12K	401001	403002/404001	401001	403002/404001	402001	402301					402601	
	ViiA7/Fast	401001	403002/404001	401001	403002/404001	402001	402301		402401	402101		402601	
	3100					402001	402301					402601	
	3130					402001	402301					402601	
3700					402001	402301					402601		
3730/3730x					402001	402301					402601		
Eppendorf®	Mastercycler®	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler® ep realplex	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler ep	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus eco	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus flat	401001		401001									
	Mastercycler nexus flat ec	401001	403002/404001	401001	403002/404001								
	Mastercycler nexus gradie	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus gradie	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus GSX1	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus GSX1	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus SX1	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler nexus SX1e	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mastercycler pro/S	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Mx3000®			401001	403002/404001	402001	402301						
	Stratagene	Mx3000P®	401001	403002/404001	401001	403002/404001	402001	402301					
		Mx3005®	401001	403002/404001	401001	403002/404001	402001	402301					
Mx4000®		401001	403002/404001	401001	403002/404001	402001	402301						
SureCycler®8800		401001	403002/404001	401001	403002/404001	402001	402301		402501	402401			
Agilent	ARIA MX G8830A										402812		
	Opticon2				403102						402101	402812	
Qiagen®	Flexigene®	401001	403002/404001	401001	403002/404001	402001	402301		402501				
	Rotor-Gene® Q	401001		401001					402501				

Brand	Instrument Model	PCR Tube		PCR Plate				
Biorad® / MJ Research®	C1000™ Touch™	401001	403002/404001	402001	402301	402501		
	DNA Engine Dyad®/Dyad	401001	403002/404001	402001	402301	402501		
	Engine Tetrad®2	401001	403002/404001	402001	402301	402501		
	MyCycler™	401001	403002/404001	402001	402301	402501		
	PTC-100®	401001	403002/404001	402001	402301	402501		
	PTC-200	401001	403002/404001	402001	402301	402501		
	PTC-225 Tetrad	401001	403002/404001	402001	402301	402501		
	S1000™	401001	403002/404001	402001	402301	402501		
	T100™	401001	403002/404001	402001	402301	402501		
	CFX Connect™	401001	403002/404001	403102	402301	402501	402101	
	CFX96™ Touch/CFX96	401001	403002/404001	403102	402301	402501	402101	402812
	Chromo4™	401001	403002/404001		402301	402501		
	DNA Engine Opticon®2	401001	403002/404001		402301	402501		
	iCycler	401001	403002/404001		402301	402501		
	iQ™5	401001	403002/404001		402301	402501		
	MiniOpticon™	401001	403002/404001		402301	402501		
	MyiQ™	401001	403002/404001		402301	402501		
	Opticon			403102		402501	402101	402812
	miniOpticon™							
	BaseStation™				402301	402501		
MJ research option							402812	
Analytik Jena / Biometra	FlexCycler	401001	403002/404001	403102	402301	402501	402101	
	T1 Thermal Cycler	401001	403002/404001	403102	402301	402501	402101	
	TGradient	401001	403002/404001	403102	402301	402501	402101	
	Tpersonal	401001	403002/404001	403102	402301	402501	402101	
	TProfessional/Basic	401001	403002/404001	403102	402301	402501	402101	
	TProfessional/Standard	401001	403002/404001	403102	402301	402501	402101	
	TProfessional/TIRO	401001	403002/404001	403102	402301	402501	402101	
	TRobot	401001	403002/404001	403102	402301	402501	402101	
	Uno	401001	403002/404001	403102	402301	402501	402101	
	Uno II	401001	403002/404001	403102	402301	402501	402101	
	q TOWER	401001	403002/404001		402301	402501		
	SpeedCycler ²	401001	403002/404001		402301	402501		
	TOptical	401001	403002/404001		402301	402501		
	Multiblock System	401001	403002/404001		402301	402501	402101	
T1 Thermal Cycler	401001	403002/404001		402301	402501	402101		
Thermo Hybaid	Omn-E	401001	403002/404001		402301	402501	402101	
	Omnigene	401001	403002/404001		402301	402501	402101	
	PCR Express	401001	403002/404001		402301	402501	402101	
	PCR Sprint	401001	403002/404001		402301	402501	402101	
	Px2	401001	403002/404001					
	Px2	401001	403002/404001		402301	402501	402101	
	Px2	401001	403002/404001		402301	402501	402101	
	Px2	401001	403002/404001		402301	402501	402101	
Touchdown	401001	403002/404001		402301	402501	402101		
3Prime	401001	403002/404001		402301	402501	402101		

Brand	Instrument Model	PCR Tube				PCR Plate			
		401001	403002/404001						
Techne®	3PrimeG	401001	403002/404001						
	3PrimeX	401001	403002/404001						
	Cyclogene™	401001	403002/404001		402001	402301	402501		
	Flexigene™	401001	403002/404001		402001	402301	402501		
	Genus	401001	403002/404001		402001	402301	402501		
	Genus Quad	401001	403002/404001		402001	402301	402501		
	Genus(TC412)	401001	403002/404001		402001	402301	402501	402101	
	Prime	401001	403002/404001		402001	402301	402501		
	Prime Elite	401001	403002/404001		402001	402301	402501		
	Prime Elite Satellite	401001	403002/404001		402001	402301	402501		
	PrimeG	401001	403002/404001		402001	402301	402501		
	Touchgene®Gradient(TC5)	401001	403002/404001	403102	402001	402301	402501	402101	
	Touchgene®X	401001	403002/404001	403102	402001	402301	402501	402101	
	PrimeQ	401001	403002/404001		402001	402301	402501		
	Quantica®	401001	403002/404001		402001	402301	402501		
	Deltacycler I	401001	403002/404001		402001	402301		402101	
	Single Block	401001	403002/404001		402001	402301		402101	
	Twin Block	401001	403002/404001		402001	402301		402101	
	Primus 96	401001	403002/404001	403102	402001	402301	402501	402101	
	peqSTAR 96	401001	403002/404001		402001	402301	402501		
SensoQuest	401001	403002/404001		402001	402301	402501			
Stratagene	401001	403002/404001		402001	402301	402501			
TaKaRa®	TP240	401001	403002/404001		402001	402301	402501		
	TP3000	401001	403002/404001	403102	402001	402301	402501	402101	
Amersham®	MegaBACE®500					402301	402501		
GE®	MegaBACE1000					402301	402501	402201	
Beckman Coul	CEQ™					402301			
Transgenomic	WAVE®System					402301	402501		
天隆	TL988-IV							402812	
Biorad® / MJ Research®	C1000™ Touch™	401001	403002/404001		402001	402301	402501		
	DNA Engine Dyad®/Dyad	401001	403002/404001		402001	402301	402501		
	Engine Tetrad®2	401001	403002/404001		402001	402301	402501		
	MyCycler™	401001	403002/404001		402001	402301	402501		
	PTC-100®	401001	403002/404001		402001	402301	402501		
	PTC-200	401001	403002/404001		402001	402301	402501		
	PTC-225 Tetrad	401001	403002/404001		402001	402301	402501		
	S1000™	401001	403002/404001		402001	402301	402501		
	T100™	401001	403002/404001		402001	402301	402501		
	CFX Connect™	401001	403002/404001	403102	402001	402301	402501	402101	
	CFX96™ Touch	401001	403002/404001	403102	402001	402301	402501	402101	
	Chromo4™	401001	403002/404001		402001	402301	402501		
	DNA Engine Opticon®2	401001	403002/404001		402001	402301	402501		
	iCycler	401001	403002/404001		402001	402301	402501		
	iQ™5	401001	403002/404001		402001	402301	402501		
	MiniOpticon™	401001	403002/404001		402001	402301	402501		


Brand	Instrument Model	PCR Tube		PCR Plate					
Analytik Jena / Biometra	MyIQ™	401001	403002/404001		402001	402301	402501		
	Opticon			403102			402501	402101	
	BaseStation™					402301	402501		
	Opticon2			403102			402501	402101	
	FlexCycler	401001	403002/404001	403102	402001	402301	402501	402101	
	T1 Thermal Cycler	401001	403002/404001	403102	402001	402301	402501	402101	
	TGradient	401001	403002/404001	403102	402001	402301	402501	402101	
	Tpersonal	401001	403002/404001	403102	402001	402301	402501	402101	
	TProfessional/Basic	401001	403002/404001	403102	402001	402301	402501	402101	
	TProfessional/Standard	401001	403002/404001	403102	402001	402301	402501	402101	
	TProfessional/TIRO	401001	403002/404001	403102	402001	402301	402501	402101	
	TRobot	401001	403002/404001	403102	402001	402301	402501	402101	
	Uno	401001	403002/404001	403102	402001	402301	402501	402101	
	Uno II	401001	403002/404001	403102	402001	402301	402501	402101	
q.TOWER	401001	403002/404001		402001	402301	402501			
SpeedCycler ²	401001	403002/404001		402001	402301	402501			
TOptical	401001	403002/404001		402001	402301	402501			
Thermo Hydraid	Multiblock System	401001	403002/404001		402001	402301	402501	402101	
	T1 Thermal Cycler	401001	403002/404001		402001	402301	402501	402101	
	Omn-E	401001	403002/404001		402001	402301	402501	402101	
	Omnigene	401001	403002/404001		402001	402301	402501	402101	
	PCR Express	401001	403002/404001		402001	402301	402501	402101	
	PCR Sprint	401001	403002/404001						
	Px2	401001	403002/404001		402001	402301	402501	402101	
	PXE	401001	403002/404001		402001	402301	402501	402101	
	Touchdown	401001	403002/404001		402001	402301	402501	402101	
	3Prime	401001	403002/404001						
	3PrimeG	401001	403002/404001						
	3PrimeX	401001	403002/404001						
	Cyclogene™	401001	403002/404001		402001	402301	402501		
	Flexigene™	401001	403002/404001		402001	402301	402501		
Techné®	Genius	401001	403002/404001		402001	402301	402501		
	Genius Quad	401001	403002/404001		402001	402301	402501		
	Genius(TC412)	401001	403002/404001		402001	402301	402501	402101	
	Prime	401001	403002/404001		402001	402301	402501		
	Prime Elite	401001	403002/404001		402001	402301	402501		
	Prime Elite Satellite	401001	403002/404001		402001	402301	402501		
	PrimeG	401001	403002/404001		402001	402301	402501		
	Touchgene®Gradient(TC)	401001	403002/404001	403102	402001	402301	402501	402101	
	Touchgene®X	401001	403002/404001	403102	402001	402301	402501	402101	
	PrimeQ	401001	403002/404001		402001	402301	402501		
	Quantica®	401001	403002/404001		402001	402301	402501		

ELISA Plates



Cat.No.	Description	Spec (well)	Color	/Pack	/Case
504201	Detachable, High Binding, Clear	96	White Frame & Transparent Well	5	50
514201	Undetachable, High Binding, Clear	96	Transparent	5	50

Cuvettes



Cat.No.	Description	/Pack	/Case
370011	3.5 mL Standard Cuvette, 10 mm, PS, Range 340-800 nm	100	1000
370211	1.5 mL Semi-Micro Cuvette, 10 mm, PS, Range 340-800 nm	100	1000



ISO 9001



ISO 11137



ISO 13485



FDA



CE

NEST[®]

Edition : First edition of March 2021

Wuxi NEST Biotechnology Co., Ltd.

JoJo Life Science UG (haftungsbeschränkt)

Biberstraße 32

89537 Giengen

Tel. 07322-9111329

Mail: info@jojo-ls.de

Web: www.jojo-ls.de

