

Power Up Your Research

- Color touch screen control panel
- Compact space-saving design
- Ideally suited for DNA and protein electrophoresis applications
- Four sets of output terminals
- Compatible with all popular gel tanks

Basic Power Supply 500V / 400mA / 120W

High Current Power Supply 300V / 2000mA / 300W

Ideal for:

- Horizontal DNA tanks
- Vertical mini protein
- Mini wet blotting



Ideal for:

- Horizontal DNA tanks
- Vertical protein tanks
- Wet and semi dry blotting
- Running multiple tanks of higher amperage



Voltage Output	500 V (1 Volt increments)
Current Output	400 mA (1 mA increments)
Output Control	Constant voltage/current
Maximum Power	120 W
Output Terminals	4 pairs
Electrical Input	115V/60Hz or 230V/50Hz
Weight/Dimensions	2.5 Lbs, 1.14 Kg / 3.5 x 5 x 8.5 inches

Voltage Output	300 V (1 Volt increments)
Current Output	2000 mA (5 mA increments)
Output Control	Constant voltage/current
Maximum Power	300 W
Output Terminals	4 pairs
Electrical Input	115V/60Hz or 230V/50Hz
Weight/Dimensions	2.7 Lbs, 1.22 Kg / 3.5 x 5 x 8.5 inches



Menge	Kurzbezeichnung	Art.-Nr.:	Listenpreis*
1 Stück	Basic Power Supply, 500 V, 400 mA, 120 W	171960	574,00 €/St.
1 Stück	High Current Power Supply, 300 V, 2000 mA, 300 W	171961	1.177,00 €/St.



PowerVolt™ Touch Screen Electrophoresis Power Supplies

JoJo
Life Science



JoJo Life Science UG (haftungsbeschränkt) - Biberstraße 32 - 89537 Giengen
Tel. 07322-9111329 - Mail: info@jojo-ls.de - Web: www.jojo-ls.de

www.jojo-ls.de – info@jojo-ls.de