

JoJo Gel-Doc-Box Basic

Photo-Gel System is a user-friendly digital system with high resolutions for gel documentation.

With a resolution up to 16 Mpixel, acquired image are of superb quality and very clear.

It is a "stand-alone" system, the use of PC is not required and all the functions are managed by the practical and powerful camera, thanks to a simply interface and the small TFT Lcd display.

Ecomic and compact system, but with powerful performance.


FRAMEWORK:

- Handle Light stainless steel pyramid darkroom with footprint 32x24cm with lateral knob to facilitate handling and displacement.
- Small front door to aid gel introduction.
- Camera adapter with filter.
- Non-slip rubber gasket under darkroom







Photo-Gel System is suitable for all type of gels and can be used in all the transilluminators available in commerce.

OPTIONAL ACCESSORIES:

- Thermal printer CP 730 for direct printing from camera.
- Gel-Quant analysis software



Photogel system Basic can be used in conjunction with our UVIFOR M20 transilluminator



JoJo Life Science UG (haftungsbeschränkt) - Biberstraße 32 - 89537 Giengen
Tel. 07322-9111329 - Mail: info@jojo-ls.de - Web: www.jojo-ls.de